


# BioBag 751: Category B 95kPa Specimen Bag



**Code: 751**

**External Dimensions**  
(once folded & sealed):  
240mm x 190mm\*

**Internal Dimensions:**  
230mm x 160mm\*

\*+/- 3mm tolerance

## Instructions:

Insert specimens and the absorbent material in the bag.

Remove the covering tape to expose the adhesive tape.

Fold at slit and firmly press.

Insert the requisition form in the back pouch.

## Barcode:

Unique barcodes for identification and traceability.

**Air Sea Containers Ltd**

521 Cavendish Street, Merseyside, CH41 8FZ, UK

▲ Tel +44 (0)151 653 1500 ▲ Email [uk@airseadg.com](mailto:uk@airseadg.com) ▲ Visit [www.airseadg.com](http://www.airseadg.com)